

## Operating manual for PY-61 series Pico second LASER

### (1) How to turn on the LASER

- ① Turn the cooling water valve.
- ② Turn on the leakage break and three-phase switch.
- ③ Turn on the main power switch.
- ④ Turn on the LASER AC switch. (**check the circulation of cooling water**)
- ⑤ Turn LASER power key. (**check lamps in securities**)
- ⑥ Push the dump reset. (**It will return to it's own position**)
- ⑦ switch the '*clock mode*' to the '*manual mode*'.
- ⑧ Push the '*Charge*' switch and then push the '*Fire*' switch.
- ⑨ Check **SIMMER** lamp. it should be turned on.
- ⑩ Repeat **step ⑧** 2 or 3 times and then change the '*manual mode*' to the '*clock mode*' after this, launch the LASER.
- ⑪ LASER shutter on after 20 minutes later.

# Operating manual for PY-61 series Pico second LASER

## (2) How to turn off

- ① Off the LASER shutter.
- ② Make the position of clock/manual switch centered.
- ③ Off the LASER key.
- ④ Off AC power.
- ⑤ Off main power switch.
- ⑥ off the leakage break.
- ⑦ off the three-phase switch.
- ⑧ turn off the water.

※ Leakage break and 3-phase switch are resided in the behind of Power supply.

